



EPOR Antibody, FITC conjugated

Product Code	CSB-PA007744LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19235
Immunogen	Recombinant Human Erythropoietin receptor protein (320-506AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Erythropoietin receptor (EPO-R), EPOR
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	EPOR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.