



EXOSC10 Antibody, Biotin conjugated

Product Code	CSB-PA007889LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01780
Immunogen	Recombinant Human Exosome component 10 protein (587-762AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exosome component 10 (EC 3.1.13.-) (Autoantigen PM/Scl 2) (P100 polymyositis-scleroderma overlap syndrome-associated autoantigen) (Polymyositis/scleroderma autoantigen 100 kDa) (PM/Scl-100) (Polymyositis/scleroderma autoantigen 2), EXOSC10, PMSCL PMSCL2 RRP6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	EXOSC10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.