



FGL1 Antibody, HRP conjugated

Product Code	CSB-PA008653LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08830
Immunogen	Recombinant Human Fibrinogen-like protein 1 protein (179-282AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibrinogen-like protein 1 (HP-041) (Hepassocin) (Hepatocyte-derived fibrinogen-related protein 1) (HFREP-1) (Liver fibrinogen-related protein 1) (LFIRE-1), FGL1, HFREP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.