



FUT3 Antibody, Biotin conjugated

Product Code	CSB-PA009077LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21217
Immunogen	Recombinant Human Galactoside 3(4)-L-fucosyltransferase protein (83-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galactoside 3(4)-L-fucosyltransferase (EC 2.4.1.65) (Blood group Lewis alpha-4-fucosyltransferase) (Lewis FT) (Fucosyltransferase 3) (Fucosyltransferase III) (FucT-III), FUT3, FT3B LE
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FUT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.