




GALR1 Antibody

Product Code	CSB-PA009223LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47211
Immunogen	Recombinant Human Galanin receptor type 1 protein (293-349AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galanin receptor type 1 (GAL1-R) (GALR-1), GALR1, GALNR GALNR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GALR1

Image



Western Blot
 Positive WB detected in: HepG2 whole cell lysate
 All lanes: GALR1 antibody at 2µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 39 kDa
 Observed band size: 39 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.