



GDF10 Antibody, HRP conjugated

Product Code	CSB-PA009343LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55107
Immunogen	Recombinant Human Growth/differentiation factor 10 protein (369-478AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth/differentiation factor 10 (GDF-10) (Bone morphogenetic protein 3B) (BMP-3B) (Bone-inducing protein) (BIP), GDF10, BMP3B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GDF10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.