



GGCX Antibody, Biotin conjugated

Product Code	CSB-PA009388LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38435
Immunogen	Recombinant Human Vitamin K-dependent gamma-carboxylase protein (383-526AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitamin K-dependent gamma-carboxylase (EC 4.1.1.90) (Gamma-glutamyl carboxylase) (Peptidyl-glutamate 4-carboxylase) (Vitamin K gamma glutamyl carboxylase), GGCX, GC
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	GGCX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.