




GPR143 Antibody

Product Code	CSB-PA009758LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51810
Immunogen	Recombinant Human G-protein coupled receptor 143 protein (314-404AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled receptor 143 (Ocular albinism type 1 protein), GPR143, OA1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPR143

Image



Western Blot

Positive WB detected in: Rat brain tissue

All lanes: GPR14 antibody at 3.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.