



HSD11B2 Antibody, Biotin conjugated

Product Code	CSB-PA010765LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80365
Immunogen	Recombinant Human Corticosteroid 11-beta-dehydrogenase isozyme 2 protein (212-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Corticosteroid 11-beta-dehydrogenase isozyme 2 (EC 1.1.1.-) (11-beta-hydroxysteroid dehydrogenase type 2) (11-DH2) (11-beta-HSD2) (11-beta-hydroxysteroid dehydrogenase type II) (11-HSD type II) (11-beta-HSD type II) (NAD-dependent 11-beta-hydroxysteroid dehydrogenase) (11-beta-HSD) (Short chain dehydrogenase/reductase family 9C member 3), HSD11B2, HSD11K SDR9C3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSD11B2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.