



HSD17B1 Antibody, Biotin conjugated

Product Code	CSB-PA010766LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14061
Immunogen	Recombinant Human Estradiol 17-beta-dehydrogenase 1 protein (268-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Estradiol 17-beta-dehydrogenase 1 (EC 1.1.1.62) (17-beta-hydroxysteroid dehydrogenase type 1) (17-beta-HSD 1) (20 alpha-hydroxysteroid dehydrogenase) (20-alpha-HSD) (E2DH) (Placental 17-beta-hydroxysteroid dehydrogenase) (Short chain dehydrogenase/reductase family 28C member 1), HSD17B1, E17KSR EDH17B1 EDH17B2 EDHB17 SDR28C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSD17B1