



HSD17B1 Antibody, Biotin conjugated

Product Code	CSB-PA010766LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14061
Immunogen	Recombinant Human Estradiol 17-beta-dehydrogenase 1 protein (268-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Estradiol 17-beta-dehydrogenase 1 (EC 1.1.1.62) (17-beta-hydroxysteroid dehydrogenase type 1) (17-beta-HSD 1) (20 alpha-hydroxysteroid dehydrogenase) (20-alpha-HSD) (E2DH) (Placental 17-beta-hydroxysteroid dehydrogenase) (Short chain dehydrogenase/reductase family 28C member 1), HSD17B1, E17KSR EDH17B1 EDH17B2 EDHB17 SDR28C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSD17B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.