



IFNGR2 Antibody, Biotin conjugated

Product Code	CSB-PA011052LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38484
Immunogen	Recombinant Human Interferon gamma receptor 2 protein (28-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon gamma receptor 2 (IFN-gamma receptor 2) (IFN-gamma-R2) (Interferon gamma receptor accessory factor 1) (AF-1) (Interferon gamma receptor beta-chain) (IFN-gamma-R-beta) (Interferon gamma transducer 1), IFNGR2, IFNGT1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNGR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.