




IRX3 Antibody

Product Code	CSB-PA011835LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78415
Immunogen	Recombinant Human Iroquois-class homeodomain protein IRX-3 protein (285-378AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iroquois-class homeodomain protein IRX-3 (Homeodomain protein IRXB1) (Iroquois homeobox protein 3), IRX3, IRXB1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	IRX3

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, Rat lung tissue

All lanes: IRX3 antibody at 5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 53 kDa

Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.