



ITGA3 Antibody, FITC conjugated

Product Code	CSB-PA011866LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26006
Immunogen	Recombinant Human Integrin alpha-3 protein (168-276AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-3 (CD49 antigen-like family member C) (FRP-2) (Galactoprotein B3) (GAPB3) (VLA-3 subunit alpha) (CD antigen CD49c) [Cleaved into: Integrin alpha-3 heavy chain; Integrin alpha-3 light chain], ITGA3, MSK18
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.