



MSN Antibody, FITC conjugated

Product Code	CSB-PA01189C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26038
Immunogen	Recombinant Human Moesin protein (466-572AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Moesin (Membrane-organizing extension spike protein), MSN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MSN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.