



KCNAB1 Antibody, FITC conjugated

Product Code	CSB-PA012013LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14722
Immunogen	Recombinant Human Voltage-gated potassium channel subunit beta-1 protein (287-401AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-gated potassium channel subunit beta-1 (EC 1.1.1.-) (K(+)) channel subunit beta-1) (Kv-beta-1), KCNAB1, KCNA1B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNAB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.