



KCNQ3 Antibody, FITC conjugated

Product Code	CSB-PA012091LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43525
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily KQT member 3 protein (763-872AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily KQT member 3 (KQT-like 3) (Potassium channel subunit alpha KvLQT3) (Voltage-gated potassium channel subunit Kv7.3), KCNQ3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNQ3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.