



GDI1 Antibody, Biotin conjugated

Product Code	CSB-PA01229D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31150
Immunogen	Recombinant Human Rab GDP dissociation inhibitor alpha protein (35-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab GDP dissociation inhibitor alpha (Rab GDI alpha) (Guanosine diphosphate dissociation inhibitor 1) (GDI-1) (Oligophrenin-2) (Protein XAP-4), GDI1, GDIL OPHN2 RABGDIA XAP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GDI1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.