



KIR2DL1 Antibody, HRP conjugated

Product Code	CSB-PA012352LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43626
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DL1 protein (162-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 2DL1 (CD158 antigen-like family member A) (MHC class I NK cell receptor) (Natural killer-associated transcript 1) (NKAT-1) (p58 natural killer cell receptor clones CL-42/47.11) (p58 NK receptor CL-42/47.11) (p58.1 MHC class-I-specific NK receptor) (CD antigen CD158a), KIR2DL1, CD158A NKAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR2DL1