



KIR2DL1 Antibody, HRP conjugated

Product Code	CSB-PA012352LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43626
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DL1 protein (162-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 2DL1 (CD158 antigen-like family member A) (MHC class I NK cell receptor) (Natural killer-associated transcript 1) (NKAT-1) (p58 natural killer cell receptor clones CL-42/47.11) (p58 NK receptor CL-42/47.11) (p58.1 MHC class-I-specific NK receptor) (CD antigen CD158a), KIR2DL1, CD158A NKAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR2DL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.