



# KLF11 Antibody, Biotin conjugated

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA012385LD01HU   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| Uniprot No.         | O14901   |
| Immunogen           | Recombinant Human Krueppel-like factor 11 protein (107-205AA)  |
| Raised In           | Rabbit   |
| Species Reactivity  | Human  |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | Biotin   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4   |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | Krueppel-like factor 11 (Transforming growth factor-beta-inducible early growth response protein 2) (TGFB-inducible early growth response protein 2) (TIEG-2), KLF11, FKLf TIEG2 |
| Immunogen Species   | Homo sapiens (Human)   |
| Research Area       | Epigenetics and Nuclear Signaling  |
| Target Names        | KLF11  |