



LTK Antibody, Biotin conjugated

Product Code	CSB-PA013230LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29376
Immunogen	Recombinant Human Leukocyte tyrosine kinase receptor protein (20-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte tyrosine kinase receptor (EC 2.7.10.1) (Protein tyrosine kinase 1), LTK, TYK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LTK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.