



LTK Antibody, Biotin conjugated

Product Code	CSB-PA013230LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29376
Immunogen	Recombinant Human Leukocyte tyrosine kinase receptor protein (20-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte tyrosine kinase receptor (EC 2.7.10.1) (Protein tyrosine kinase 1), LTK, TYK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LTK