



MAGEA3 Antibody, HRP conjugated

Product Code	CSB-PA013329LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43357
Immunogen	Recombinant Human Melanoma-associated antigen 3 protein (2-104AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-associated antigen 3 (Antigen MZ2-D) (Cancer/testis antigen 1.3) (CT1.3) (MAGE-3 antigen), MAGEA3, MAGE3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MAGEA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.