



MAP3K4 Antibody, HRP conjugated

Product Code	CSB-PA013426LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6R4
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 4 protein (53-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 4 (EC 2.7.11.25) (MAP three kinase 1) (MAPK/ERK kinase kinase 4) (MEK kinase 4) (MEKK 4), MAP3K4, KIAA0213 MAPKKK4 MEKK4 MTK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.