



MAP3K4 Antibody, Biotin conjugated

Product Code	CSB-PA013426LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6R4
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 4 protein (53-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 4 (EC 2.7.11.25) (MAP three kinase 1) (MAPK/ERK kinase kinase 4) (MEK kinase 4) (MEKK 4), MAP3K4, KIAA0213 MAPKKK4 MEKK4 MTK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K4