



MAP3K6 Antibody, Biotin conjugated

Product Code	CSB-PA013428LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95382
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 6 protein (905-990AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 6 (EC 2.7.11.25) (Apoptosis signal-regulating kinase 2), MAP3K6, ASK2 MAPKKK6 MEKK6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAP3K6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.