



MASP2 Antibody, HRP conjugated

Product Code	CSB-PA013510LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00187
Immunogen	Recombinant Human Mannan-binding lectin serine protease 2 protein (189-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mannan-binding lectin serine protease 2 (EC 3.4.21.104) (MBL-associated serine protease 2) (Mannose-binding protein-associated serine protease 2) (MASP-2) [Cleaved into: Mannan-binding lectin serine protease 2 A chain; Mannan-binding lectin serine protease 2 B chain], MASP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MASP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.