



MC1R Antibody, Biotin conjugated

Product Code	CSB-PA013558LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01726
Immunogen	Recombinant Human Melanocyte-stimulating hormone receptor protein (1-37AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanocyte-stimulating hormone receptor (MSH-R) (Melanocortin receptor 1) (MC1-R), MC1R, MSHR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MC1R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.