



MC1R Antibody, Biotin conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA013558LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q01726 |
| Immunogen | Recombinant Human Melanocyte-stimulating hormone receptor protein (1-37AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Melanocyte-stimulating hormone receptor (MSH-R) (Melanocortin receptor 1) (MC1-R), MC1R, MSHR |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | MC1R |