



MCHR2 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA013585LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q969V1 |
| Immunogen | Recombinant Human Melanin-concentrating hormone receptor 2 protein (1-39AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Melanin-concentrating hormone receptor 2 (MCH receptor 2) (MCH-R2) (MCHR-2) (G-protein coupled receptor 145) (GPRv17) (MCH-2R) (MCH2) (MCH2R), MCHR2, GPR145 SLT |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | MCHR2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |