



METAP2 Antibody, HRP conjugated

Product Code	CSB-PA013716LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50579
Immunogen	Recombinant Human Methionine aminopeptidase 2 protein (313-449AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methionine aminopeptidase 2 (MAP 2) (MetAP 2) (EC 3.4.11.18) (Initiation factor 2-associated 67 kDa glycoprotein) (p67) (p67eIF2) (Peptidase M), METAP2, MNPEP P67EIF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	METAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.