



MGST1 Antibody, FITC conjugated

Product Code	CSB-PA013791LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10620
Immunogen	Recombinant Human Microsomal glutathione S-transferase 1 protein (63-96AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microsomal glutathione S-transferase 1 (Microsomal GST-1) (EC 2.5.1.18) (Microsomal GST-I), MGST1, GST12 MGST
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MGST1