



MST1 Antibody, HRP conjugated

Product Code	CSB-PA015055LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26927
Immunogen	Recombinant Human Hepatocyte growth factor-like protein (363-465AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatocyte growth factor-like protein (Macrophage stimulatory protein) (Macrophage-stimulating protein) (MSP) [Cleaved into: Hepatocyte growth factor-like protein alpha chain; Hepatocyte growth factor-like protein beta chain], MST1, D3F15S2 DNF15S2 HGFL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MST1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.