



MSX1 Antibody, Biotin conjugated

Product Code	CSB-PA015068LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28360
Immunogen	Recombinant Human Homeobox protein MSX-1 protein (19-179AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein MSX-1 (Homeobox protein Hox-7) (Msh homeobox 1-like protein), MSX1, HOX7
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MSX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.