



MSX2 Antibody, FITC conjugated

Product Code	CSB-PA015069LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35548
Immunogen	Recombinant Human Homeobox protein MSX-2 protein (50-120AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein MSX-2 (Homeobox protein Hox-8), MSX2, HOX8
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MSX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.