



MYH7 Antibody, HRP conjugated

Product Code	CSB-PA015300LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12883
Immunogen	Recombinant Human Myosin-7 protein (1802-1915AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-7 (Myosin heavy chain 7) (Myosin heavy chain slow isoform) (MyHC-slow) (Myosin heavy chain, cardiac muscle beta isoform) (MyHC-beta), MYH7, MYHCB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.