



MYT1 Antibody, Biotin conjugated

Product Code	CSB-PA015387LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01538
Immunogen	Recombinant Human Myelin transcription factor 1 protein (317-451AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin transcription factor 1 (MyT1) (Myelin transcription factor I) (MyTI) (PLPB1) (Proteolipid protein-binding protein), MYT1, KIAA0835 KIAA1050 MTF1 MYTI PLPB1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MYT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.