



NMT1 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA015900LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P30419 |
| Immunogen | Recombinant Human Glycylpeptide N-tetradecanoyltransferase 1 protein (1-158AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Glycylpeptide N-tetradecanoyltransferase 1 (EC 2.3.1.97) (Myristoyl-CoA:protein N-myristoyltransferase 1) (NMT 1) (Type I N-myristoyltransferase) (Peptide N-myristoyltransferase 1), NMT1, NMT |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cancer |
| Target Names | NMT1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |