



NOVA1 Antibody, HRP conjugated

Product Code	CSB-PA015957LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51513
Immunogen	Recombinant Human RNA-binding protein Nova-1 protein (6-174AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein Nova-1 (Neuro-oncological ventral antigen 1) (Onconeural ventral antigen 1) (Paraneoplastic Ri antigen) (Ventral neuron-specific protein 1), NOVA1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NOVA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.