



OCRL Antibody, Biotin conjugated

Product Code	CSB-PA016267LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01968
Immunogen	Recombinant Human Inositol polyphosphate 5-phosphatase OCRL-1 protein (121-271AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol polyphosphate 5-phosphatase OCRL-1 (EC 3.1.3.36) (Lowe oculocerebrorenal syndrome protein), OCRL, INPP5F OCRL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OCRL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.