



OCRL Antibody, Biotin conjugated

Product Code	CSB-PA016267LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01968
Immunogen	Recombinant Human Inositol polyphosphate 5-phosphatase OCRL-1 protein (121-271AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol polyphosphate 5-phosphatase OCRL-1 (EC 3.1.3.36) (Lowe oculocerebrorenal syndrome protein), OCRL, INPP5F OCRL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OCRL