



OPRM1 Antibody, HRP conjugated

Product Code	CSB-PA016361LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35372
Immunogen	Recombinant Human Mu-type opioid receptor protein (1-68AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mu-type opioid receptor (M-OR-1) (MOR-1) (Mu opiate receptor) (Mu opioid receptor) (MOP) (hMOP), OPRM1, MOR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OPRM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.