



PCSK6 Antibody, HRP conjugated

Product Code	CSB-PA017645LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29122
Immunogen	Recombinant Human Proprotein convertase subtilisin/kexin type 6 protein (596-710AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proprotein convertase subtilisin/kexin type 6 (EC 3.4.21.-) (Paired basic amino acid cleaving enzyme 4) (Subtilisin-like proprotein convertase 4) (SPC4) (Subtilisin/kexin-like protease PACE4), PCSK6, PACE4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PCSK6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.