



PEX10 Antibody, FITC conjugated

Product Code	CSB-PA017794LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60683
Immunogen	Recombinant Human Peroxisome biogenesis factor 10 protein (63-157AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisome biogenesis factor 10 (Peroxin-10) (Peroxisomal biogenesis factor 10) (Peroxisome assembly protein 10) (RING finger protein 69), PEX10, RNF69
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PEX10