



PEX5 Antibody, Biotin conjugated

Product Code	CSB-PA017805LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50542
Immunogen	Recombinant Human Peroxisomal targeting signal 1 receptor protein (149-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal targeting signal 1 receptor (PTS1 receptor) (PTS1R) (PTS1-BP) (Peroxin-5) (Peroxisomal C-terminal targeting signal import receptor) (Peroxisome receptor 1), PEX5, PXR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PEX5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.