



PF4 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA017809LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P02776 |
| Immunogen | Recombinant Human Platelet factor 4 protein (32-101AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Platelet factor 4 (PF-4) (C-X-C motif chemokine 4) (Iroplact) (Oncostatin-A) [Cleaved into: Platelet factor 4, short form (Endothelial cell growth inhibitor)], PF4, CXCL4 SCYB4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | PF4 |