



PF4 Antibody, Biotin conjugated

Product Code	CSB-PA017809LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02776
Immunogen	Recombinant Human Platelet factor 4 protein (32-101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet factor 4 (PF-4) (C-X-C motif chemokine 4) (Iroplact) (Oncostatin-A) [Cleaved into: Platelet factor 4, short form (Endothelial cell growth inhibitor)], PF4, CXCL4 SCYB4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PF4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.