



PNLIP Antibody, HRP conjugated

Product Code	CSB-PA018262LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16233
Immunogen	Recombinant Human Pancreatic triacylglycerol lipase protein (60-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pancreatic triacylglycerol lipase (PL) (PTL) (Pancreatic lipase) (EC 3.1.1.3), PNLIP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PNLIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.