



POLD1 Antibody, FITC conjugated

Product Code	CSB-PA018303LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28340
Immunogen	Recombinant Human DNA polymerase delta catalytic subunit protein (1-123AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase delta catalytic subunit (EC 2.7.7.7) (EC 3.1.11.-) (DNA polymerase subunit delta p125), POLD1, POLD
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.